

GALAXY S

Electrosurgical Generator

U_{max} 0,5 kVp max 60 Watt 5 CH. O101 101 10 USES L CABLE = 3 M L₁ = 360 MM

| | | | |
|-----------|--|-----------|--|
| 77-074-51 | | 77-076-51 | |
| 77-075-51 | | 77-077-51 | |

66-403-00 GALAXY SMART UTERIS UNIT

66-403--08 PURE SALINE UTERIS BUNDLE



67-901-12 GALAXY SMART FOOT SWITCH

67-591-45 CONNECTING CABLE FOR ALL PLASMA SALINE RESECTOSCOPES

SYSTEM XS SYSTEM S SYSTEM M SYSTEM L



GYNIX™ is a Registered Trademark of Schultz Medical Ltd.

SCHULTZ MEDICAL™
THE SINGLE-USE INSTRUMENT COMPANY

Schultz Medical Ltd, Unit 6,
Commerce Park, Crowland Street,
Southport PR9 7RL
+44 (0)1772 250 298
info@schultzmedical.co.uk
www.schultzmedical.co.uk

Europe Office:
Alt-Heerdt 104,
40549 Dusseldorf,
Germany
info@schultz-medical.com

USA Office:
45 Rockefeller Plaza,
20th Floor, New York, NY 10111
Tel +1 212 332 3215
info@schultz-medical.com



Comfort, Care & Trust

GALAXY S Electrosurgical Generator

SCHULTZ MEDICAL™
THE SINGLE-USE INSTRUMENT COMPANY

GALAXY S **NEW**

Electrosurgical Generator

Special Electrosurgical Unit For Pure Bipolar Applications:

- Hystero-Resection
- Bipolar Needle Electrodes




SALINE
0,9% NaCl




PURE PLASMA SALINE GENERATOR

SCHULTZ MEDICAL
THE SINGLE-USE INSTRUMENT COMPANY



EASY AND INTUITIVE TOUCH SCREEN HANDLING



REAL TIME OUTPUT POWER MONITORING



AUTOMATIC POWER REGULATION
RESULT-ORIENTED POWER MANAGEMENT



ADDITIONAL INFORMATION ON THE CONNECTED INSTRUMENT SAVED ON THE DIGITAL CHIP



QUICK RELEASE BUTTON FOR RECOMMENDED STANDARD SETTINGS



SDS CONNECTION WITH AUTOMATIC INSTRUMENT RECOGNITION
THE SDS CABLES COMMUNICATE WITH THE GALAXY UNIT IN AN INTERACTIVE WAY